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
In all exercises, the emphasis should be on stretching. Following the exercises, remember to maintain good posture, i.e. to sit with the breastbone raised.

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Mobilization of the spinal column, circulation enhancement
Rock pelvis back and forth

Mobilization of the spinal column, circulation enhance
Move pelvis from side to side

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
Stabilization of the spinal column
1. Alternately: Raise one leg and opposing arm diagonally
2. Alternately: Raise one leg, turn trunk in the direction of the raised leg; maintain good posture!

Mobility of upper spinal column with arm movement
Bend and stretch trunk with folded arms, emphasis on stretching.

Mobilization, stability, stretching
1. soft backwards motion of alternating arms, eyes follow arm movement (maintain good posture!)
2. Standing, arms rotated to the outside, palms downward on Wedge Ball Cushion, lower arm stretching.

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Relaxation, circulation, venous gymnastics, stimulation of foot reflex zones - without shoes - alternately: press heel downwards and pull toes upwards, then press toes downwards and pull heels upwards.

Find balance when seated
Balance, back away from support, feet suspended

Balance, "waking up" the body at intervals
Stand on the Wedge Ball Cushion with no shoes; at first with support, then try to balance without hand support

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Filling

Deflate

Fold over the half of the **AIRGO-Active Cushion** with the valve (1), insert the lightly moistened tip of the needle into the valve and allow air to escape by pinching it gently (it is not necessary to screw the pump onto the needle for this purpose).

Inflate

Attach the needle (stored in the pump handle) to the end of the pump and moisten the tip of the needle slightly. Fold over the half of the **AIRGO-Active Cushion** with the valve (1), insert the tip of the needle into the valve and pump (2).

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
3) Just right/Gonflage optimum

Test for correct inflation volume:

Press with your hand on one half of the **AIRGO-Active Cushion**. The inflation volume is just right when the other half remains flat (i.e. does not bulge or go slack).

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Position I / Placement I

Correct
Bon

Tip the seat of your chair backwards slightly. Place the **AIRGO Active Cushion** about 5 cm away from the back of the chair on the seat. Make sure your knees are slightly above your buttocks when you are seated. Avoid tipping the seat surface towards the front (see figure at the right).

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
Position II

Place the **AIRGO-Active Cushion** about 5 cm away from the back of the chair on its seat (valve pointing downwards).

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Care - Use a mild agent for cleaning plastic surfaces such as Eco-Power-Cleaner for gentle cleaning of the **AIRGO-Active Cushion**. Avoid contact with sharp or pointed objects. Use only the pump supplied with the cushion together with its needle attachment to avoid damage to the valve.

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The **AIRGO-Active Cushion** assumes body temperature after only a few minutes. The odour of the new **AIRGO-Active Cushion** dissipates within a short time as in a new car.

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GB Important information

The Dynair® Premium is shipped from our factory with a ready-to-use, comfortable air filling. You can, however, use the enclosed TOGU® air pump for easy regulation of the air filling.

The inflation / deflation needle is recessed into the head of the air pump handle under the small cap. To **deflate** the Dynair® Premium somewhat, moisten the needle, then introduce it straight into the valve (see Fig.) and press out the excess air.

To **inflate** the Dynair® Premium somewhat, screw the needle onto the air pump and moisten it, then introduce it straight into the valve (see Fig.) and pump the needed air into it.

Important: Before introducing the needle pinch the Dynair® Premium valve between two fingers and pull up a little (see Fig.), since otherwise you might puncture the other side of the product. To avoid damaging the valve, use only the pump supplied with the product.

Never pump up the Dynair® Premium until it is filled taut. Less internal air pressure is more beneficial, since only then can the desired fine vibrations be transmitted to the vertebral disc and column.

Product care: Use a mild agent for cleaning plastics, e.g. Eco-Power Cleaner, for gentle cleaning of the Dynair® Premium.

Avoid contact with sharp or pointed objects.

The Dynair® Premium assumes body temperature within a few minutes. The nobbly side provides for enhanced air circulation while sitting.

Sewing Instructions for Ball Cushion Cover:

Fabric required: For Dynair® Premium Kids cover approx. 55 cm in diameter. For Dynair® Premium cover approx. 61 cm in diameter

Draw a circle in one of the following sizes: For Dynair® Premium Kids: 52 cms in diameter (radius 26

cms). For Dynair® Premium: 60 cms in diameter (radius 30 cms)

The simplest way to draw a circle: Fix a thread or string (radius see above) to a pencil. Keep hold of the end of the string with a finger and draw the circle with the other hand. This is how you obtain the pattern.

Fold fabric at marked line, fold in approx. 0,5 cm and stitch around narrow edge.

Pull in elastic band loosely. Place fabric around cushion. Tighten elastic band until fabric fits tightly around cushion without any creases. Shorten length of elastic band to fit and sew two ends of elastic band together. Close drawstring seam.

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