Anleitung

2607203

Pedalo Stabilisator „Profi“

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14003000

Pedalo® Stabilizer 'Professional'

Before using the device for the first time
(or before starting with assembly),
please **read instruction completely** and keep it

Holz-Hoerz GmbH
Dottinger Str. 71
D-Münsingen

Tel.: (+49) 07381 / 9357-0
Fax (+49) 07381 / 9357-40

eMail: info@pedalo.de
<http://www.pedalo.de>

Important information and warning hints:

Warning! This device is a very demanding one! It is to be used with caution, since it requires great skill, so as to avoid falls or collisions causing injury to the user or third parties. Especially when used for the first times and/or when doing exercises in standing position, a special degree in caution, attention, circumspection and a high degree in motor abilities is assumed. In case of doubt first use the device only with assistance and/or appropriate protective equipment.

On request place the Pedalo-Stabilizer in front of a wall bar or the like.
When used by children the supervision by an adult is recommended.

Warning! The packing is not part of the product and have to be disposed.

Keep plastic bags and films out of the reach of babies and young children – **RISK OF SUFFOCATION!**

Warning! The device is only to be used for the purposes described in this user manual.

Warning! Only use device on level ground with even and firm surfaces. On sensible grounds use an anti-slip mat under it if necessary. The Pedalo Stabilizer may only be used by one person!

Protect item from damp and wet.

Maintenance and cleaning: Can be wiped off with a damp cloth. Do not use abrasive cleaning agents.

The maximum weight capacity is 150 kg with static and even loading of all 4 springs. In case of dynamic or one-sided loading the value is significantly reduced.

This device is meant for using from persons at the age of 6 up.

Check the screw, plug and mount connections before every use. Also check the springs and metal parts every four weeks for deformations.

It is never permissible to use the Stabilizer when it shows visible damages.

Depending on using intensity the springs have to be replaced every two years at least. Please always replace all springs completely. Defective and worn parts should be replaced by original replacement parts immediately.

Assembly of the frame:

Put together the two footrests (with the black plastic feet) and the two upper brackets (with the pre-assembled eyebolts).

Now a safety bolt is inserted through the holes for each pipe connection.


Note that the plastic part is securely engaged with the pipe.

Assembly of the ropes and springs: The 4 springs are prepared as shown

1. The orange safety rope is attached as shown with a rope loop on the small screw carabiner. The clamp of the safety rope should be as close as possible to the carabiner.

2. Now push the safety rope through the inside of the spring. The folding carabiner is hung down into the safety rope.

Attention! The closed spring eye must be below.

3. Now the springs are attached to the stainless steel suspension bar. The suspension bar is inserted through the opening of the carabiner. The spring is attached to the first notch of the suspension bar.

4. Then hang both suspension bars with the springs on the hooks in the tube frame. Pull all safety ropes down.

Variant 1 - without spring function

Now a standing platform is attached (middle notch at the suspension angle). Make sure that the rubber o-rings are below the carabiners. They prevent the position of the carabiner from being inadvertently dislodged.

With this construction the platforms will swing horizontally in all directions while doing exercises.

Variant 2 - with spring function


The platform is now hung on the springs and safety ropes. Now the platform moves horizontally and vertically during exercises.

The safety rope prevents the springs from being overloaded and the floorboard can touch the floor.

If both small standing platforms are to be connected to a large standing platform, proceed as follows:

First of all, place all 4 springs in the upper suspension completely outwards. Place the springs in the corner of the small standing platforms in the middle catch. Now insert the studs of the split standing platform of the metal brackets of the standing platform. In the middle, the standing platform is only laid together like a puzzle piece. Now you can put the Pedalo® Anti-Slip-Mat on the platform.

The supplied Anti-Slip-Mat is perfect for a safe stand on the device. The mat can be wiped off with a damp cloth. Please note that the mat should always be on the Stabilizer only for the duration of the use. For storage and cleaning please remove the mat from the standing platform.

Change of vibration:

By changing the suspension of the spring, quite different application requirements can be produced. The springs can be adjusted above or below by one or more catches.

Make sure that all spring pairs are at the top and bottom in the same pegs.

The standing platform could otherwise hang obliquely or possibly individual springs could be loaded too heavily.

Attention! Before you get on, please check whether all the springs are correctly installed in the notches. If necessary, move the standing platform slightly back and forth with your hand.

Attention!
After each change, the requirements for the users can change significantly. Therefore, mount the device carefully. Possibly with support or suitable holding device such as a rung wall.

Change of tipping behavior:

Decentralized suspension of the standing platforms for strengthening the outer or inner tapes and correction of the leg axis.

Even with exercises on the individual standing platforms it is possible to adjust the spring suspension completely according to your wishes.

If a small standing platform is loaded with full body weight, it is necessary to hang this on 4 springs.

Either with the two springs of the second small standing platform or with additional springs including carabiner.

Item 627042: Tension spring stainless steel

Item 627108: Carabiner

Accessory: Pedalo®-Footboards Seesaw Item 653101

Increases the level of difficulty.

Aims on foot and leg axis stability.

Pedalo Footboards seesaw offer a linear movement left-right. Can be attached to other Pedalo products but are also suitable to be used individually on the floor (retractable pins on the bottom side).

Accessory: Pedalo®-Softpads Item 144056

Enforce the training effect.

Specialty foam, 2 pcs. 30 x 14 x 2.5 cm